

Distinct detection of intact and defective proviral DNA for HIV subtype B and C


N.V.E.J. Buchholtz¹, M.M. Nühn¹, T.C.M. de Jong¹, A.M.J. Wensing^{1, 2, 3}, J. Symons¹, M. Nijhuis^{1, 2, 3}

1. Virology, Department of Medical Microbiology, University Medical Center Utrecht (UMCU), Utrecht, the Netherlands
 2. Wits Reproductive Health and HIV Institute (Wits RHI), University of the Witwatersrand, Johannesburg, South Africa
 3. Ndlovu Research Consortium, Elandsdoorn, South Africa

Introduction


- The HIV reservoir is the main obstacle to a cure
- During treatment the reservoir can be replication competent (intact, 2.4%) or replication incompetent (defective, 97.6%)
- Reservoir quantification is essential to monitor cure interventions

To achieve a global cure it is important to quantify the reservoir of the most dominant HIV subtypes, B and C



Methods

- Multiplex Digital Droplet PCR
- 2 primer/probe sets used to quantify the reservoir



Optimized protocol to distinguish intact and defective proviral HIV subtype B and C

- 1 fluorescent signal → defective proviral DNA (*Psi* or *Env*)
- 2 fluorescent signals → intact proviral DNA



Results

- Equal efficacy of the subtype B and C primers and probes was defined on Gblocks (fig.1)
- Limit of detection, which defines the sensitivity of the assay was determined at 7 copies with >95% certainty (fig.2)
- Limit of blank, which defines the specificity was determined at 0.90 copies at 60°C and 1.23 copies at 55°C (fig.2&3)
- Natural variation of 1 nucleotide change in *Env* can still be detected by lowering the annealing temperature to 55°C (fig.3)




Fig 1. Detection efficacy of the *psi* and *env* region with both subtype B and C primers/probes on Gblocks.

copy input	% of positive signals	mean copy output	standard deviation
6	97,5	6.38	3,68
5	92,5	4.22	2,97
4	90,0	3.50	2,36
3	95,1	3.21	2,22
2	83,3	1.91	1,63
1	61,4	1.11	1,19
0	14,0	0.12	0,48

LoB = mean_{blank} + 1.645 * (SD blank) LoD = LoB + 1.645 * (SD sample 95% positivity)
 LoB = 0.12 + 1.645 * 0.48 = 0.90 LoD = 0.90 + (1.645 * 3.68) = 6.95

Fig 2 Specificity and sensitivity of the subtype B and C primers and probes at 60°C. The Limit of Blank (LoB) was determined at 0.90 copies and the limit of detection (LoD) at 7 copies.

Copies/22µl	Annealing temperature			
	60°C	58 °C	55 °C	53 °C
Env intact	1293,6	1238,6	1001,0	1273,1
Env probe G-A mutation 5th nucleotide	2,0	2,8	1021,4	1372,1
Env probe G-A mutation 13th nucleotide	0,6	0,0	840,4	1177,7
Env hypermutation probe	0,0	0,9	0,4	2,3

LoB = 0.18 + 1.645 * 0.64 = 1.23

Fig 3. Detection of the subtype B and C primers/probes at the *Env* region with lower annealing temperatures, tested on Gblocks. The specificity of the assay at 55°C was determined at 1.23 copies.

Discussion

- This IPDA is a specific and sensitive assay for the detection of subtype B and C defective and intact proviral DNA
- Decreasing the annealing temperature to 55°C enables correct quantification of patient isolates with 1 nucleotide change in the *Env* gene
- This IPDA is most likely also able to correctly quantify subtype D and F, since these strains have no sequence differences in *psi* and *env*

The IPDA can form the basis for large-scale evaluation of global cure interventions in clinical studies